
Laboratory
Equipment for ICSI



Vardhman Medicare Pvt. Ltd.

ICSI System



*RI Integra 3 Micro-manipulator
with Nikon Inverted Research
Microscope (Ti2-U)*

1. Integrated with:

- Make RI, model Integra 3 Micromanipulator;
- Nikon Inverted Research Microscope (model Ti2U);
- Large built-in mechanical XY stage
- 3.5'' color LCD touchscreen display for temperature control

2. Thermosafe™ Air Heating System with Heated Metal Insert: 25mm aperture (standard); 16mm aperture (optional); ITO Heated Glass Insert (optional);

3. Motion sensor LED light underneath the stage for improved visibility when changing objectives;

4. New extra-fine mechanisms for smooth and precise movement of the stage and micropipettes.



Reach out to us at:

L-90/164, 1st and 2ng Floor, Street number 7,
Mahipalpur Extension, New Delhi, 110037,
India.

Contact: 9818245106 / 9811226146

Email ID: sales@vardhmanmedicare.net

Website: <http://vardhmanmedicare.net/>



Vardhman Medicare Pvt. Ltd.